






The next generation of male sling technology™ *only from AMS*



Innovative Procedure

- New transobturator approach for treating mild to moderate SUI
- The sling effects a proximal relocation of the urethral bulb which causes a circumferential coaptation of the sphincteric urethra
- The force of the sling is applied parallel to the urethral lumen, unlike compressive devices which apply force into the lumen

Promising Results

- Including the feasibility study, nearly 200 men have been treated with AdVance™, with as much as 36 months of follow up data
- AdVance has been successful in treating mild to moderate SUI with no retention or erosion

Learn more at www.AMSAdVance.com



The next generation of male sling technology™ *only from AMS*

New mesh and needle technology!

All AdVance™ components have been designed specifically for use with the male anatomy



Patient Selection

- Mild to moderate SUI
- 1-5 pads per day
- Post RP or TURP
- Residual sphincter function
- Prior radiation is ok

Procedural Information

- Minimally invasive
- Outpatient procedure
- 20-30 minute case time
- Spinal or general anesthesia
- Does not preclude future AUS or penile prosthesis

Learn more at www.AMSAdVance.com



American Medical Systems, Inc.
World Headquarters
10700 Bren Road West
Minnetonka, MN 55343
USA
Phone: 952 930 6000
Fax: 952 930 6157
www.AmericanMedicalSystems.com

For a complete list of indications, contraindications, warnings and precautions, refer to the product Instructions for Use or contact your AMS representative.